



Featureset Reference Manual

SVCam-U3V Series

Digital USB3Vision Area Scan Cameras

Version 1.7.2 Date 2020/11/09

exo

Device Control

Category for Device information and control.

Feature

Device Vendor Name
Type: String
 GenICam Name: DeviceVendorName
 Name of the manufacturer of the device.

Feature

Device Model Name
Type: String
 GenICam Name: DeviceModelName
 Model name of the device.

Feature

Device Manufacturer Info
Type: String
 GenICam Name: DeviceManufacturerInfo
 Manufacturer information about the device.

Feature

Device Version
Type: String
 GenICam Name: DeviceVersion
 Version of the device.

Feature

DeviceFirmwareVersion
Type: String
 GenICam Name: DeviceFirmwareVersion
 Device Firmware Version

Feature

Device ID
Type: String
 GenICam Name: DeviceID
 Device identifier (serial number).

Feature

Device User ID
Type: String
 GenICam Name: DeviceUserID
 User-programmable device identifier.

Feature

Device Temperature Selector
Type: Enumeration
 GenICam Name: DeviceTemperatureSelector
 Selects the location within the device, where the temperature will be measured.

Enumeration Entities

Name	GenICam Name
Mainboard	Mainboard
Description:	This enumeration value selects the temperature measured on the mainboard.

Feature

Device Temperature

Type: Float

GenICam Name: DeviceTemperature

Device temperature in degrees Celsius (C). It is measured at the location selected by DeviceTemperatureSelector.

Feature

Device Reset

Type: Command

GenICam Name: DeviceReset

Resets the device to its power up state.

Feature

U3V Device Capability

Type: Integer

GenICam Name: U3VDeviceCapability

Describes Device Capability.

Feature

User Name Available

Type: Boolean

GenICam Name: UserNameAvailable

Set if User Defined Name is available.

Feature

Access Privilege Available

Type: Boolean

GenICam Name: AccessPrivilegeAvailable

Set if Heartbeat/Access Privilege is available.

Feature

Message Channel Supported

Type: Boolean

GenICam Name: MessageChannelSupported

Set if the device supports a Message Channel.

Feature

Timestamp Supported

Type: Boolean

GenICam Name: TimestampSupported

Set if the device supports a timestamp register.

Feature

String Encoding

Type: Integer

GenICam Name: StringEncoding

String Encoding of the BRM.

Feature

Family Register Available

Type: Boolean

GenICam Name: FamilyRegisterAvailable

Set if the device supports the Family Name register.

Feature

SBRM Supported

Type: Boolean

GenICam Name: SBRMSupported

Set if the device supports a SBRM.

Feature

Endianess Registers SupportedType: Boolean

GenICam Name: EndianessRegistersSupported

Set if the device supports the Protocol Endianess and Implementation Endianess registers.

Feature

Written Length Field SupportedType: Boolean

GenICam Name: WrittenLengthFieldSupported

Set to 1 if the device sends the length_written field in the SCD section of the WriteMemAck command.

Image Format Control

This category includes items that control the size of the acquired image and the format of the transferred pixel data

Feature

Sensor Width Type: Integer

GenICam Name: SensorWidth
This is a read only element. It is an integer that indicates the actual width of the camera

Feature

Sensor Height Type: Integer

GenICam Name: SensorHeight
This is a read only element. It is an integer that indicates the actual height of the camera

Feature

X Offset Type: Integer

GenICam Name: OffsetX
This value sets the left offset for the area of interest in pixels, i.e., the distance in pixels between the left side of the sensor and the left side of the image.

Feature

Y Offset Type: Integer

GenICam Name: OffsetY
This value sets the top offset for the area of interest, i.e., the distance in pixels between the top of the sensor and the top of the image,

Feature

Width Type: Integer

GenICam Name: Width
This value sets the width of the area of interest in pixels.

Feature

Height Type: Integer

GenICam Name: Height
This value sets the height of the area of interest in pixels.

Feature

Max Width Type: Integer

GenICam Name: WidthMax
This is a read only element. It is an integer that indicates maximum allowed width of the image in pixels taking into account any function that may limit the allowed width.

Feature

Max Height Type: Integer

GenICam Name: HeightMax
This is a read only element. It is an integer that indicates maximum allowed height of the image in pixels taking into account any function that may limit the allowed height.

Feature

Pixel Size

Type: Enumeration

GenICam Name: PixelSize

This is a read only feature. This enumeration provides a list of values that indicate the depth of the pixel values in the acquired images in bits per pixel.

Enumeration Entities

Name	GenICam Name
Bpp8	Bpp8
Description:	<i>This enumeration value indicates that the depth of the pixel values in the acquired images is 8 bits per pixel.</i>
Bpp12	Bpp12
Description:	<i>This enumeration value indicates that the depth of the pixel values in the acquired images is 12 bits per pixel.</i>
Bpp16	Bpp16
Description:	<i>This enumeration value indicates that the depth of the pixel values in the acquired images is 16 bits per pixel.</i>

Feature

Pixel Format

Type: Enumeration

GenICam Name: PixelFormat

Format of the pixel to use for acquisition. It represents all the informations provided by PixelCoding, PixelSize, PixelColorFilter but combined in one single value.

Enumeration Entities

Name	GenICam Name
Mono8	Mono8
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Mono 8'.</i>
Mono12Packed	Mono12Packed
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Mono 12 packed'.</i>
Mono16	Mono16
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Mono 16'.</i>
BayerGR8	BayerGR8
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GR 8'.</i>
BayerRG8	BayerRG8
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer RG 8'.</i>
BayerGB8	BayerGB8
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GB 8'.</i>
BayerBG8	BayerBG8
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer BG 8'.</i>
BayerGR12Packed	BayerGR12Packed
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GR 12 packed'.</i>

Enumeration Entities

Name	GenICam Name
BayerRG12Packed	BayerRG12Packed
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer RG 12 packed'.</i>
BayerGB12Packed	BayerGB12Packed
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GB 12 packed'.</i>
BayerBG12Packed	BayerBG12Packed
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer BG 12 packed'.</i>
BayerGR12	BayerGR12
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GR 12'.</i>
BayerRG12	BayerRG12
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer RG 12'.</i>
BayerGB12	BayerGB12
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer GB 12'.</i>
BayerBG12	BayerBG12
Description:	<i>This enumeration value sets the pixel format of the image data transmitted by the camera to 'Bayer BG 12'.</i>

Feature

Pixel Color Filter

Type: Enumeration

GenICam Name: PixelColorFilter

This is a read only feature. This enumeration provides a list of values that indicate the alignment of the camera's Bayer filter to the pixels in the acquired images.

Enumeration Entities

Name	GenICam Name
None	None
Description:	<i>This enumeration value indicates that no Bayer filter is present on the camera.</i>
BayerRG	BayerRG
Description:	<i>This enumeration value indicates that the Bayer filter has an RG/GB alignment to the pixels in the acquired images.</i>
BayerGR	BayerGR
Description:	<i>This enumeration value indicates that the Bayer filter has an RG/GB alignment to the pixels in the acquired images.</i>
BayerGB	BayerGB
Description:	<i>This enumeration value indicates that the Bayer filter has an RG/GB alignment to the pixels in the acquired images.</i>
BayerBG	BayerBG
Description:	<i>This enumeration value indicates that the Bayer filter has an RG/GB alignment to the pixels in the acquired images.</i>



Feature

PixelDynamicRangeMin

Type: Integer

GenICam Name: PixelDynamicRangeMin
Minimum value that can be returned during the digitization process. This corresponds to the darkest value of the camera.

Feature

PixelDynamicRangeMax

Type: Integer

GenICam Name: PixelDynamicRangeMax
Maximum value that will be returned during the digitization process. This corresponds to the brightest value of the camera.

Feature

Binning Horizontal

Type: Integer

GenICam Name: BinningHorizontal
This enumeration controls the horizontal binning setting

Feature

Binning Vertical

Type: Integer

GenICam Name: BinningVertical
This enumeration controls the horizontal binning setting

Feature

Decimation Horizontal

Type: Enumeration

GenICam Name: DecimationHorizontal
This enumeration reduces the horizontal resolution(Width) of the image by the specified horizontal decimation factor.No increased frame rate

Enumeration Entities

Name	GenICam Name
Off	Off
On	On

Feature

Decimation Vertical

Type: Enumeration

GenICam Name: DecimationVertical
This enumeration reduces the vertical resolution(Height) of the image by the specified vertical decimation factor.Increased frame rate

Enumeration Entities

Name	GenICam Name
Off	Off
On	On

Feature

ReverseX

Type: Boolean

GenICam Name: ReverseX
Flip horizontally the image sent by the device. The Region of interest is applied after the flipping.

Feature

ReverseY

Type: Boolean

GenlCam Name: ReverseY

Flip vertically the image sent by the device. The Region of interest is applied after the flipping.

Acquisition Control

This category includes items used to set the image acquisition parameters and to start and stop acquisition

Feature

Acquisition Mode

Type: Enumeration

GenICam Name: AcquisitionMode
This enumeration sets the image acquisition mode.

Enumeration Entities

Name	GenICam Name
Single Frame	SingleFrame
Description:	This enumeration value sets the camera's acquisition mode to single frame
Multi Frame	MultiFrame
Description:	This enumeration value sets the camera's acquisition mode to multi frame.
Continuous	Continuous
Description:	This enumeration value sets the camera's acquisition mode to continuous.

Feature

Acquisition Start

Type: Command

GenICam Name: AcquisitionStart
Starts the Acquisition of the device. The number of frames captured is specified by AcquisitionMode.

Feature

Acquisition Stop

Type: Command

GenICam Name: AcquisitionStop
Stops the Acquisition of the device at the end of the current Frame. It is mainly used when AcquisitionMode is Continuous but can be used in any acquisition mode.

Feature

Trigger Selector

Type: Enumeration

GenICam Name: TriggerSelector
This enumeration lists the types of trigger that are available for selection. Once a trigger type has been selected, all the other trigger features will be applied to the selected trigger.

Enumeration Entities

Name	GenICam Name
Acquisition Start	AcquisitionStart
Description:	This enumeration value selects the acquisition start trigger.
Frame Start	FrameStart
Description:	This enumeration value selects the frame start trigger.

Feature

Trigger Mode

Type: Enumeration

GenICam Name: TriggerMode
This enumeration provides a list of the values available for setting the trigger mode for the selected trigger.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value sets the mode for the selected trigger to off.</i>
On	On
Description:	<i>This enumeration value sets the mode for the selected trigger to on.</i>

Feature

Generate Software Trigger

Type: Command

GenICam Name: TriggerSoftware

This command generates a software trigger signal. The software trigger signal will be used if the trigger source is set to 'software'.

Feature

Trigger Source

Type: Enumeration

GenICam Name: TriggerSource

This enumeration lists the available trigger sources for the selected trigger.

Enumeration Entities

Name	GenICam Name
Trigger Software	Software
Description:	<i>This enumeration value sets the source for the selected trigger to software trigger.</i>
Line 1	Line1
Description:	<i>This enumeration value sets the source for the selected trigger to line 1.</i>

Feature

Trigger Activation

Type: Enumeration

GenICam Name: TriggerActivation

This enumeration lists the trigger activation types available for the selected trigger.

Enumeration Entities

Name	GenICam Name
Rising Edge	RisingEdge
Description:	<i>This enumeration value sets the trigger to be valid when the trigger signal is going high.</i>
Falling Edge	FallingEdge
Description:	<i>This enumeration value sets the trigger to be valid when the trigger signal is going low.</i>

Feature

Trigger Delay

Type: Float

GenICam Name: TriggerDelay

Trigger Delay specifies the absolute delay in microseconds (us) to apply the trigger reception before effectively activating it.

Feature

Exposure Mode

Type: Enumeration

GenICam Name: ExposureMode

This enumeration lists the available exposure modes.



Enumeration Entities

Name	GenICam Name
Timed	<i>Timed</i>
Description:	<i>This enumeration value sets the exposure mode to timed.</i>
Trigger Width	<i>TriggerWidth</i>
Description:	<i>This enumeration value sets the exposure mode to trigger width.</i>

Feature

Acquisition Frame Rate

Type: Float

GenICam Name: AcquisitionFrameRate

This float value sets the camera's acquisition frame rate in Hz.

Feature

Exposure Time

Type: Float

GenICam Name: ExposureTime

This float value sets the camera's exposure time in microseconds.

Feature

Exposure Auto

Type: Enumeration

GenICam Name: ExposureAuto

Sets the automatic exposure mode when ExposureMode is Timed. The exact algorithm used to implement this control is device-specific.

Enumeration Entities

Name	GenICam Name
Off	<i>Off</i>
Description:	<i>Exposure duration is user controlled using ExposureTime.</i>
Once	<i>Once</i>
Description:	<i>Exposure duration is adapted once by the device. Once it has converged, it returns to the Off state.</i>
Continuous	<i>Continuous</i>
Description:	<i>Exposure duration is constantly adapted by the device to maximize the dynamic range.</i>

Feature

Exposure First

Type: Boolean

GenICam Name: ExposureAutoOrder

This value sets the priority of both exposure and gain settings. True means exposure before gain.

Feature

ExposureTimeMin

Type: Float

GenICam Name: ExposureTimeMin

This float value sets the camera's minimum exposure time for auto exposure in microseconds.

Feature

ExposureTimeMax

Type: Float

GenICam Name: ExposureTimeMax

This float value sets the camera's maximum exposure time for auto exposure in microseconds.



Feature

Sensor Shutter Mode Type: Enumeration

GenICam Name: SensorShutterMode
This enumeration sets the shutter mode.

Enumeration Entities

Name	GenICam Name
Global Shutter	Global
Description:	<i>This enumeration value sets the camera's shutter mode to global shutter</i>
Rolling Shutter	Rolling
Description:	<i>This enumeration value sets the camera's shutter mode to rolling shutter.</i>
Global Reset Shutter	GlobalReset
Description:	<i>This enumeration value sets the camera's shutter mode to global reset shutter.</i>

Feature

PIV Mode Type: Enumeration

GenICam Name: PIVMode
This enumeration activates the PIV mode.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration disables the PIV mode</i>
On	On
Description:	<i>This enumeration enables the PIV mode</i>

Analog Control

This category includes items that control the analog characteristics of the video signal

Feature

Gain Selector

Type: Enumeration

GenICam Name: GainSelector

This enumeration selects the gain control to configure. Once a gain control has been selected, all changes to the gain settings will be applied to the selected control.

Enumeration Entities

Name	GenICam Name
All	All
Description:	<i>This enumeration value selects all available gain controls for adjustment.</i>

Feature

Gain (dB)

Type: Float

GenICam Name: Gain

Sets the dB value of the selected gain control.

Feature

Digital Gain

Type: Integer

GenICam Name: DigitalGain

Sets the dB value of the selected digital gain control.

Feature

Black Level Selector

Type: Enumeration

GenICam Name: BlackLevelSelector

This enumeration selects the black level control to configure. Once a black level control has been selected, all changes to the black level settings will be applied to the selected control.

Enumeration Entities

Name	GenICam Name
All	All
Description:	<i>This enumeration value selects all available black level controls for adjustment.</i>

Feature

Offset(Raw)

Type: Integer

GenICam Name: BlackLevelRaw

This value sets the selected black level control as an integer.

Feature

Digital Offset

Type: Integer

GenICam Name: DigitalOffset

This value sets the control of the digital Offset.

Feature

Gain Auto

Type: Enumeration

GenICam Name: GainAuto

Sets the automatic gain control (AGC) mode. The exact algorithm used to implement AGC is device-specific.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	Gain is User controlled using Gain.
Once	Once
Description:	Gain is automatically adjusted once by the device. Once it has converged, it automatically returns to the Off state.
Continuous	Continuous
Description:	Gain is constantly adjusted by the device.

Feature

Gain Speed Type: Enumeration
GenICam Name: GainSpeed
This value sets the speed of the autogain:standard speed or fast speed.

Enumeration Entities

Name	GenICam Name
Standard	Standard
Fast	Fast

Feature

GainAutoMin Type: Float
GenICam Name: GainAutoMin
This Float value sets the minimum applied gain for AGC in dB

Feature

GainAutoMax Type: Float
GenICam Name: GainAutoMax
This Float value sets the maximum applied gain for AGC in dB

Feature

Autogain Level Type: Integer
GenICam Name: GainAutoLevel
this feature sets the automatic gain level

Feature

Balance Ratio Selector Type: Enumeration
GenICam Name: BalanceRatioSelector
This enumeration selects a balance ratio control to configuration. Once a balance ratio control has been selected, all changes to the balance ratio settings will be applied to the selected control.

Enumeration Entities

Name	GenICam Name
Red	Red
Description:	This enumeration value selects the red balance ratio control for adjustment.
Green	Green
Description:	This enumeration value selects the green balance ratio control for adjustment.

Enumeration Entities

Name	GenICam Name
Blue	Blue
Description:	<i>This enumeration value selects the blue balance ratio control for adjustment.</i>

Feature

Balance Ratio

Type: Float

GenICam Name: BalanceRatio

Controls ratio of the selected color component to a reference color component. It is used for white balancing.

Feature

Gain Auto Balance

Type: Enumeration

GenICam Name: GainAutoBalance

Sets the mode for automatic gain balancing between the aps. The gain coefficients of each channel or tap are adjusted so they are matched.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>Gain tap balancing is user controlled using Gain.</i>
Once	Once
Description:	<i>Gain tap balancing is automatically adjusted once by the device. Once it has converged, it automatically returns to the Off state.</i>
Continuous	Continuous
Description:	<i>Gain tap balancing is constantly adjusted by the device.</i>
Reset	Reset
Description:	<i>Gain tap balancing is resetted.</i>

Feature

BalanceWhiteAuto

Type: Enumeration

GenICam Name: BalanceWhiteAuto

Controls the mode for automatic white balancing between the color channels.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>Auto white balancing is user controlled using Gain.</i>
Once	Once
Description:	<i>Auto white balancing is automatically adjusted once by the device. Once it has converged, it automatically returns to the Off state.</i>
Continuous	Continuous
Description:	<i>Auto white balancing is constantly adjusted by the device.</i>
Reset	Reset
Description:	<i>Auto white balancing is resetted.</i>

Digital IO Control

This category includes items used to control the operation of the camera's digital I/O lines

Feature

LineSelector

Type: Enumeration

GenICam Name: LineSelector

This enumeration selects the I/O line to configure.

Enumeration Entities

Name	GenICam Name
Output 0 - Line0	Line0
Description:	<i>This enumeration value selects 'Output 0' for configuration.</i>
Output 1 - Line1	Line1
Description:	<i>This enumeration value selects 'Output 1' for configuration.</i>
Output 2 - Line2	Line2
Description:	<i>This enumeration value selects 'Output 2' for configuration.</i>
Output 3 - Line3	Line3
Description:	<i>This enumeration value selects 'Output 3' for configuration.</i>
Output 4 - Line4	Line4
Description:	<i>This enumeration value selects 'Output 4' for configuration.</i>
Uart In - Line5	Line5
Description:	<i>This enumeration value selects 'Uart In' for configuration.</i>
Trigger - Line6	Line6
Description:	<i>This enumeration value selects 'Trigger' for configuration.</i>
Sequencer - Line7	Line7
Description:	<i>This enumeration value selects 'Sequencer' for configuration.</i>
Debouncer - Line8	Line8
Description:	<i>This enumeration value selects 'Debouncer' for configuration.</i>
Prescaler - Line9	Line9
Description:	<i>This enumeration value selects 'Prescaler' for configuration.</i>
Logic A - Line15	Line15
Description:	<i>This enumeration value selects 'Logic A' for configuration.</i>
Logic B - Line16	Line16
Description:	<i>This enumeration value selects 'Logic B' for configuration.</i>
Lens TXD - Line17	Line17
Description:	<i>This enumeration value selects 'Lens TXD' for configuration.</i>
Pulse 0 - Line18	Line18
Description:	<i>This enumeration value selects 'Pulse 0' for configuration.</i>
Pulse 1 - Line19	Line19
Description:	<i>This enumeration value selects 'Pulse 1' for configuration.</i>
Pulse 2 - Line20	Line20
Description:	<i>This enumeration value selects 'Pulse 2' for configuration.</i>



Enumeration Entities

Name	GenICam Name
Pulse 3 - Line21	Line21
Description:	<i>This enumeration value selects 'Pulse 3' for configuration.</i>
Uart2 In - Line22	Line22
Description:	<i>This enumeration value selects 'Uart2' In for configuration.</i>
Input 0 - Line10	Line10
Description:	<i>This enumeration value selects 'Input 0' for configuration.</i>
Input 1 - Line11	Line11
Description:	<i>This enumeration value selects 'Input 1' for configuration.</i>
Input 2 - Line12	Line12
Description:	<i>This enumeration value selects 'Input 2' for configuration.</i>
Input 3 - Line13	Line13
Description:	<i>This enumeration value selects 'Input 3' for configuration.</i>
Input 4 - Line14	Line14
Description:	<i>This enumeration value selects 'Input 4' for configuration.</i>

Feature

LineMode

Type: Enumeration

GenICam Name: LineMode

Controls if the physical Line is used to Input or Output a signal.

Enumeration Entities

Name	GenICam Name
Input	Input
Description:	<i>This enumeration value sets the mode for the selected line to 'input', i.e., the line is used to input an electrical signal.</i>
Output	Output
Description:	<i>This enumeration value sets the mode for the selected line to 'Output', i.e., the line is used to output an electrical signal.</i>

Feature

LineInverter

Type: Boolean

GenICam Name: LineInverter

Controls the inversion of the signal of the selected input or output Line.

Feature

LineStatus

Type: Boolean

GenICam Name: LineStatus

This boolean value indicates the current logical state for the selected line.

Feature

LineSource

Type: Enumeration

GenICam Name: LineSource

Selects which internal acquisition or I/O source signal to output on the selected Line.LineMode must be Output.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value sets the source signal for the selected output line to Off.</i>
Input 0	Input0
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Input 0'.</i>
Input 1	Input1
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Input 1'.</i>
Input 2	Input2
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Input 2'.</i>
Input 3	Input3
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Input 3'.</i>
Input 4	Input4
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Input 4'.</i>
User Output 0	UserOutput0
Description:	<i>This enumeration value sets the source signal for the selected output line to 'User Output 0'.</i>
User Output 1	UserOutput1
Description:	<i>This enumeration value sets the source signal for the selected output line to 'User Output 1'.</i>
User Output 2	UserOutput2
Description:	<i>This enumeration value sets the source signal for the selected output line to 'User Output 2'.</i>
User Output 3	UserOutput3
Description:	<i>This enumeration value sets the source signal for the selected output line to 'User Output 3'.</i>
User Output 4	UserOutput4
Description:	<i>This enumeration value sets the source signal for the selected output line to 'User Output 4'.</i>
UART Out	UartOut
Description:	<i>This enumeration value sets the source signal for the selected output line to 'UART Out'.</i>
Strobe 0	Strobe0
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Strobe 0'.</i>
Strobe 1	Strobe1
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Strobe 1'.</i>
Strobe 2	Strobe2
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Strobe 2'.</i>
Strobe 3	Strobe3
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Strobe 3'.</i>
PWM 0	PWM0
Description:	<i>This enumeration value sets the source signal for the selected output line to 'PWM 0'.</i>
PWM 1	PWM1
Description:	<i>This enumeration value sets the source signal for the selected output line to 'PWM 1'.</i>
PWM 2	PWM2
Description:	<i>This enumeration value sets the source signal for the selected output line to 'PWM 2'.</i>
PWM 3	PWM3
Description:	<i>This enumeration value sets the source signal for the selected output line to 'PWM 3'.</i>

Enumeration Entities

Name	GenICam Name
Expose	Expose
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Expose'.</i>
Readout	Readout
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Readout'.</i>
SeqPulse A	SeqPulseA
Description:	<i>This enumeration value sets the source signal for the selected output line to 'SeqPulse A'.</i>
SeqPulse B	SeqPulseB
Description:	<i>This enumeration value sets the source signal for the selected output line to 'SeqPulse B'.</i>
Sequencer active	SeqActive
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Sequencer active'.</i>
Debouncer	Debouncer
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Debouncer'.</i>
Prescaler	Prescaler
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Prescaler'.</i>
Logic	Logic
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Logic'.</i>
Lens Uart out	LensUartOut
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Lens Uart out'.</i>
Triggerfeedback	Triggerfeedback
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Triggerfeedback'.</i>
Uart 2 Out	Uart2Out
Description:	<i>This enumeration value sets the source signal for the selected output line to 'Uart 2 Out'.</i>

Feature

LineStatusAll

Type: Integer

GenICam Name: LineStatusAll

This integer value is a single bitfield that indicates the current logical state of all available lines at time of polling.

Feature

UserOutputSelector

Type: Enumeration

GenICam Name: UserOutputSelector

Selects which bit of the User Output register will be set by 'User Output Value'.

Enumeration Entities

Name	GenICam Name
UserOutput 0	UserOutput0
Description:	<i>This enumeration value selects user settable output 'UserOutput 0' for configuration.</i>
UserOutput 1	UserOutput1
Description:	<i>This enumeration value selects user settable output 'UserOutput 1' for configuration.</i>
UserOutput 2	UserOutput2
Description:	<i>This enumeration value selects user settable output 'UserOutput 2' for configuration.</i>

Enumeration Entities

Name	GenICam Name
UserOutput 3	UserOutput3
Description:	<i>This enumeration value selects user settable output 'UserOutput 3' for configuration.</i>
UserOutput 4	UserOutput4
Description:	<i>This enumeration value selects user settable output 'UserOutput 4' for configuration.</i>

Feature

UserOutputValue Type: Boolean
GenICam Name: UserOutputValue
This boolean value sets the state of the selected user settable output signal.

Feature

UserOutputValueAll Type: Integer
GenICam Name: UserOutputValueAll
Sets the value of all the bits of the User Output register.

Feature

UserOutputValueAllMask Type: Integer
GenICam Name: UserOutputValueAllMask
Sets the write mask to apply to the value specified by UserOutputValueAll before writing it in the User Output register.setting the user Output register using UserOutputValueAll will only change the bits that have a corresponding bit in the mask set to one.

Strobe Control

This category includes items used to set the parameters for the integrated strobe controller

Feature

Strobe Selector Type: Enumeration

GenICam Name: StrobeSelector
Selects which strobe to configure.

Enumeration Entities

Name	GenICam Name
Strobe 0	Strobe0
Description:	Selects 'Strobe 0' to configure.
Strobe 1	Strobe1
Description:	Selects 'Strobe 1' to configure.
Strobe 2	Strobe2
Description:	Selects 'Strobe 2' to configure.
Strobe 3	Strobe3
Description:	Selects 'Strobe 3' to configure.

Feature

Strobe Polarity Type: Enumeration

GenICam Name: StrobePolarity
This Enumeration sets the camera's strobe polarity.

Enumeration Entities

Name	GenICam Name
positive	positive
Description:	Sets the strobe polarity to be positive.
negative	negative
Description:	Sets the strobe polarity to be negative.

Feature

Strobe Duration Type: Float

GenICam Name: StrobeDuration
This float value sets the camera's strobe duration in microseconds.

Feature

Strobe Delay Type: Float

GenICam Name: StrobeDelay
This float value sets the camera's strobe delay in microseconds.

Enhanced IO

This category includes items used to control the integrated PWM Controller

Feature		
PWMEnable	Type: Boolean	
GenICam Name: PWMEnable		
This feature enables or disables the PWM settings.		
Feature		
PWMMax	Type: Integer	
GenICam Name: PWMMax		
This feature sets the common PWM frequency in system clock ticks.		
Feature		
PWMChange0	Type: Integer	
GenICam Name: PWMChange0		
This feature sets the Duty Cycle for PWM register		
Feature		
PWMChange1	Type: Integer	
GenICam Name: PWMChange1		
This feature sets the Duty Cycle for PWM register		
Feature		
PWMChange2	Type: Integer	
GenICam Name: PWMChange2		
This feature sets the Duty Cycle for PWM register		
Feature		
PWMChange3	Type: Integer	
GenICam Name: PWMChange3		
This feature sets the Duty Cycle for PWM register		
Feature		
SeqTrigger	Type: Command	
GenICam Name: SeqTrigger		
This feature starts the sequencer with a software signal.		
Feature		
SeqTriggermode	Type: Enumeration	
GenICam Name: SeqTriggermode		
This feature selects the sequencer trigger mode.		
<u>Enumeration Entities</u>		
Name	GenICam Name	
Trigger on high level	LevelHigh	
Description:	This feature sets the sequencer trigger on high level.	
Trigger on rising edge	RisingEdge	
Description:	This feature sets the sequencer trigger on rising edge.	

Feature

SeqSelector

GenICam Name: SeqSelector
Index to the sequencer array.

Type: Integer

Feature

SeqCount

GenICam Name: SeqCount
Number of sequencer inputs.

Type: Integer

Feature

SeqEnable

GenICam Name: SeqEnable
Activate the sequencer.

Type: Boolean

Feature

SeqLoop

GenICam Name: SeqLoop
Sequencer loop enables/disables.

Type: Boolean

Feature

DebounceDuration

GenICam Name: DebounceDuration
Sets the value for the debouncer duration.

Type: Integer

Feature

PrescaleDivisor

GenICam Name: PrescaleDivisor
Sets the value of the prescaler divisor.

Type: Integer

Feature

LogicFunction

GenICam Name: LogicFunction
Controls if the physical Line is used to Input or Output a signal.

Type: Enumeration

Enumeration Entities

Name	GenICam Name
AND	AND_Function
Description: AND logic block.	
OR	OR_Function
Description: OR logic block.	
XOR	XOR_Function
Description: XOR logic block.	
NAND	NAND_Function
Description: NAND logic block.	
NOR	NOR_Function
Description: NOR logic block.	

Enumeration Entities

Name	GenICam Name
XNOR	<i>XNOR_Function</i>
Description:	<i>XNOR logic block.</i>
Trigger Enable	<i>TRIGGER_ENABLE</i>
Description:	<i>Trigger enable logic block.</i>

Feature

SeqInterval	Type: Integer
GenICam Name: SeqInterval	
Sets the duration of the sequencer segments.	

Feature

SeqPulseAStart	Type: Integer
GenICam Name: SeqPulseAStart	
Begin of the integration time in the sequencer segment.	

Feature

SeqPulseAStop	Type: Integer
GenICam Name: SeqPulseAStop	
End of the integration time in the sequencer segment.	

Feature

SeqPulseBStart	Type: Integer
GenICam Name: SeqPulseBStart	
Begin of the PWM mask in the sequencer segment.	

Feature

SeqPulseBStop	Type: Integer
GenICam Name: SeqPulseBStop	
End of the integration time in the sequencer segment.	

LUT Control

Category that includes the LUT control features.

Feature

LUT Selector

Type: Enumeration

GenICam Name: LUTSelector
Selects which LUT to control.

Enumeration Entities

Name	GenICam Name
Luminance	Luminance
Description:	Selects the Luminance LUT.

Feature

LUT Enable

Type: Boolean

GenICam Name: LUTEnable
This boolean value enables the selected LUT.

Feature

LUT Index

Type: Integer

GenICam Name: LUTIndex
Control the index (offset) of the coefficient to access in the selected LUT.

Feature

LUT Value

Type: Integer

GenICam Name: LUTValue
Returns the Value at entry LUTIndex of the LUT selected by LUTSelector.

Feature

Gamma

Type: Float

GenICam Name: Gamma
Controls the gamma correction of pixel intensity. This is typically used to compensate for nonlinearity of the display system (such as CRT). The realization of the gamma correction is implemented using a LUT. Therefore, some LUT functionality is not available when gamma correction is activated.

Customer ID Protection

This category includes items used for licensing.

Feature

Customer ID

Type: Integer

GenICam Name: CustomerID

This integer value includes the customer id set by manufacturer.

Feature

Customer Value

Type: Integer

GenICam Name: CustomerValue

This integer value includes the customer value set by customer.

Feature

Customer Value Key

Type: Integer

GenICam Name: CustomerValueKey

This integer value includes the key for setting the customer value by customer.

Feature

Customer Data Index

Type: Integer

GenICam Name: CustomerDataIndex

Control the index (offset) of the CustomerData Array.

Feature

Customer Data Value

Type: Integer

GenICam Name: CustomerDataValue

Returns the Value at entry CustomerDataIndex of the CustomerData Array.

Lens Control

This category includes items used to set the parameters for the MFT Lens Control

Feature

Lens Set State Type: Enumeration

GenICam Name: MFTLensState
This enumeration activate/deactivate the MFT lens.

Enumeration Entities

Name	GenICam Name
inactive	inactive
Description:	This enumeration value sets the lens inactive.
active	active
Description:	This enumeration value activates the lens.

Feature

Lens Name Type: String

GenICam Name: MFTLensName
Shows the MFT lens name

Feature

Focal Length Type: Integer

GenICam Name: MFTFocalLength
This integer value sets the focal length.

Feature

Focus Type: Integer

GenICam Name: MFTFocus
This integer value sets the focus.

Feature

Focus Unit Type: Enumeration

GenICam Name: MFTFocusUnit
This integer value sets the focus unit: 1mm or 1/10mm.

Enumeration Entities

Name	GenICam Name
1mm	mft_unit_mm
Description:	focus unit 1mm.
1/10mm	mft_unit_deci_mm
Description:	focus unit 1/10mm.

Feature

Aperture Type: Integer

GenICam Name: MFTAperture
This integer value sets the aperture.

Feature

Lens Control Type

GenICam Name: LensControlType
Selected Lenscontroller type.

Type: Enumeration

Enumeration Entities

Name	GenICam Name
none	none
Description:	No lens controller available.
Birger Mount	birger
Description:	Birger Mount connected to camera.
Varioptic	varioptic
Description:	Varioptic connected to camera.
SVCam-EF lens adapter	svcamef
Description:	EF lens connected to camera.

Feature

Focus

GenICam Name: LensControlFocus
This integer value sets the focus.

Type: Integer

Feature

Iris

GenICam Name: LensControlIris
This integer value sets the iris.

Type: Integer

Feature

Lens Init

GenICam Name: LensControlLensInit
Resets the Lens.

Type: Command

Defect Pixel Correction

This category includes items used to set the parameters for the Defect Pixel Correction

Feature

Control Type: Enumeration

GenICam Name: DefectPixelCorrectionEnable_Control
This enumeration provides a list of the values available for controlling the defect pixel correction.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value disables the defect pixel correction.</i>
On	On
Description:	<i>This enumeration value enables the defect pixel correction.</i>

Feature

Mark defect pixels Type: Enumeration

GenICam Name: DefectPixelCorrectionMark_Control
This enumeration provides a list of the values available for setting test modes for the defect pixel correction.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value disables the marking of the defect pixels.</i>
Mark	Mark
Description:	<i>This enumeration value enables the marking of the defect pixels.</i>

Feature

Defect Pixel Map Selector Type: Enumeration

GenICam Name: DefectPixelCorrection_MapSelect
This enumeration provides a list of maps available for the defect pixel correction.

Enumeration Entities

Name	GenICam Name
Factory Map	factory
Description:	<i>This enumeration value selects the factory map, containing defect pixels from sensor data sheet.</i>
Custom1 Map	custom1
Description:	<i>This enumeration value selects the custom1 map, handled by customer.</i>
Custom2 Map	custom2
Description:	<i>This enumeration value selects the custom2 map, handled by customer.</i>

Feature

Defect Pixel Map Max Size Type: Integer

GenICam Name: DefectPixelCorrection_MapMaxSize
This integer value reads the maximal number of defect pixels per map.

Feature

Defect Pixel Map Size Type: Integer
GenICam Name: DefectPixelCorrection_MapSize
This integer value reads the current number of defect pixels of the selected map.

Feature

Dynamic Hot Pixel Correction Control Type: Enumeration
GenICam Name: DynamicHotPixelCorrectionEnable_Control
This enumeration provides a list of the values available for controlling the dynamic hot pixel correction.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value disables the dynamic hot pixel correction.</i>
On	On
Description:	<i>This enumeration value enables the dynamic hot pixel correction.</i>

Shading

This category includes items used to set the parameters for the Shading Correction

Feature

Shading Control

Type: Enumeration

GenICam Name: Shading_Control

This enumeration provides a list of the values available for controlling the shading correction.

Enumeration Entities

Name	GenICam Name
Off	Off
Description:	<i>This enumeration value disables the shading correction.</i>
On	On
Description:	<i>This enumeration value enables the shading correction.</i>

Feature

Shading Map Selector

Type: Enumeration

GenICam Name: Shading_MapSelect

This enumeration provides a list of maps available for the shading correction.

Enumeration Entities

Name	GenICam Name
Shading Map 0	ShadingMap0
Description:	<i>This enumeration value selects the shading map 0.</i>
Shading Map 1	ShadingMap1
Description:	<i>This enumeration value selects the shading map 1.</i>
Shading Map 2	ShadingMap2
Description:	<i>This enumeration value selects the shading map 2.</i>

Transport Layer Control

Category that contains the Transport layer control features.

Feature		
Payload Size	Type: Integer	
GenICam Name: PayloadSize Provides the number of bytes transferred for each image or chunk on the stream channel. This includes any end-of-line, end-of-frame statistics or other stamp data. This is the total size of data payload for a data block.		
Feature		
GenCP Version Major	Type: Integer	
GenICam Name: GenCPVersionMajor Major version of the specification.		
Feature		
GenCP Version Minor	Type: Integer	
GenICam Name: GenCPVersionMinor Minor version of the specification.		
Feature		
Max Response Time	Type: Integer	
GenICam Name: U3VMaxDeviceResponseTime Max Resonse Time in ms.		
Feature		
U3V Version Major	Type: Integer	
GenICam Name: U3VVersionMajor U3V Version.		
Feature		
U3V Version Minor	Type: Integer	
GenICam Name: U3VVersionMinor U3V Version.		
Feature		
U3V Capability	Type: Integer	
GenICam Name: U3VCPCapability Indicates additional features on the control channel.		
Feature		
U3V SIRM Available	Type: Boolean	
GenICam Name: U3VCPSIRMAvailable Set if the device supports at least one device streaming interface.		
Feature		
U3V EIRM Available	Type: Boolean	
GenICam Name: U3VCPEIRMAvailable Set if the device supports at least one device event interface.		

Feature

U3V IIDC2 Available

GenICam Name: U3VCPIIDC2Available

Set if the device supports IIDC2 register map.

Type: Boolean

Feature

Max Command Transfer Length

GenICam Name: U3VMaxCommandTransferLength

Specifies the max suuported command transfer length of the device.

Type: Integer

Feature

Max Ack Transfer Length

GenICam Name: U3VMaxAcknowledgeTransferLength

Specifies the max suuported ack transfer length of the device.

Type: Integer

Feature

Number of Stream Channels

GenICam Name: U3VNumberOfStreamChannels

Number of Stream Channels and its Corresponding Streaming Interface Register Maps.

Type: Integer

Feature

Current Speed

GenICam Name: U3VCurrentSpeed

Specifies the current speed of the USB link.

Type: Integer

Feature

DeviceTapGeometry

GenICam Name: DeviceTapGeometry

This device tap geometry feature describes the geometrical properties characterizing the taps of a camera as presented at the output of the device.

Type: Enumeration

Enumeration Entities

Name	GenICam Name
Geometry_1X_1Y	Geometry_1X_1Y
Description: 1X 1Y tap geometry.	
Geometry_2XE_1Y	Geometry_2XE_1Y
Description: 2XE 1Y tap geometry.	
Geometry_1X_2YE	Geometry_1X_2YE
Description: 1X 2YE tap geometry.	
Geometry_2XE_2YE	Geometry_2XE_2YE
Description: 2XE 2YE tap geometry.	
Geometry_2X_1Y	Geometry_2X_1Y
Description: 2X 1Y tap geometry.	

User Set Control

User Sets provides the features used to save camera settings to on-board non-volatile memory.

Feature

UserSetSelector

Type: Enumeration

GenICam Name: UserSetSelector

This enumeration selects the user set to load, save, or configure. Once a user set has been selected, all changes to the user set settings will be applied to the selected user set.

Enumeration Entities

Name	GenICam Name
Default User Set	Default
Description:	<i>This enumeration value selects the default user set. This is a user set that contains factory settings. It is read-only and cannot be modified.</i>
User Set 1	UserSet1
Description:	<i>This enumeration value selects the user set 1 configuration set.</i>
User Set 2	UserSet2
Description:	<i>This enumeration value selects the user set 2 configuration set.</i>
User Set 3	UserSet3
Description:	<i>This enumeration value selects the user set 3 configuration set.</i>
User Set 4	UserSet4
Description:	<i>This enumeration value selects the user set 4 configuration set.</i>
User Set 5	UserSet5
Description:	<i>This enumeration value selects the user set 5 configuration set.</i>
User Set 6	UserSet6
Description:	<i>This enumeration value selects the user set 6 configuration set.</i>
User Set 7	UserSet7
Description:	<i>This enumeration value selects the user set 7 configuration set.</i>
User Set 8	UserSet8
Description:	<i>This enumeration value selects the user set 8 configuration set.</i>

Feature

UserSetLoad

Type: Command

GenICam Name: UserSetLoad

This command loads the User Set specified by UserSetSelector to the device and makes it active.

Feature

UserSetSave

Type: Command

GenICam Name: UserSetSave

This command copies the parameters in the current active user set into the selected user set in the camera's non-volatile memory.

Feature

UserSetDefault

Type: Enumeration

GenICam Name: UserSetDefault

This enumeration sets the user set to be used as the default startup set. The user set that has been selected as the default startup set will be loaded as the active set whenever the camera is powered on or reset.

Enumeration Entities

Name	GenICam Name
Default User Set	<i>Default</i>
Description:	<i>This enumeration value sets the default user set as the default startup set.</i>
User Set 1	<i>UserSet1</i>
Description:	<i>This enumeration value sets user set 1 as the default startup set.</i>
User Set 2	<i>UserSet2</i>
Description:	<i>This enumeration value sets user set 2 as the default startup set.</i>
User Set 3	<i>UserSet3</i>
Description:	<i>This enumeration value sets user set 3 as the default startup set.</i>
User Set 4	<i>UserSet4</i>
Description:	<i>This enumeration value sets user set 3 as the default startup set.</i>
User Set 5	<i>UserSet5</i>
Description:	<i>This enumeration value sets user set 5 as the default startup set.</i>
User Set 6	<i>UserSet6</i>
Description:	<i>This enumeration value sets user set 6 as the default startup set.</i>
User Set 7	<i>UserSet7</i>
Description:	<i>This enumeration value sets user set 7 as the default startup set.</i>
User Set 8	<i>UserSet8</i>
Description:	<i>This enumeration value sets user set 8 as the default startup set.</i>

Debug

Feature

RegisterAddress

Type: Integer

GenICam Name: RegisterAddress

This feature allows the direct access to the camera via register for debug purpose.

Feature

RegisterValue(DEC)

Type: Integer

GenICam Name: RegisterValue_dec

This feature outputs the value of the accessed register in decimal.

Feature

RegisterValue(HEX)

Type: Integer

GenICam Name: RegisterValue_hex

This feature outputs the value of the accessed register in hexadecimal.

Feature

Refresh Register

Type: Command

GenICam Name: Refresh

This feature allows to update the register content manually.

Feature

FB Statistic

Type: Integer

GenICam Name: FB_Statistic

Camera framebuffer informations.

Feature

FB Used Max

Type: Integer

GenICam Name: FB_UsedDW_Max

Camera framebuffer informations.

Feature

FB Discard Count

Type: Integer

GenICam Name: FB_Discard_Count

Camera framebuffer informations.

Feature

FB Critical Set

Type: Integer

GenICam Name: FB_Critical_Set

Camera framebuffer informations.

Feature

FB Critical Reset

Type: Integer

GenICam Name: FB_Critical_Reset

Camera framebuffer informations.



Feature

FB Critical Watermark	Type: Integer
GenlCam Name: FB_Critical_Watermark	
Camera framebuffer informations.	